MINIASPEED Battery Pro

GUARANTEED

3

YEARS

HOSPITAL

ISTRUZIONI D'USO

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INSTRUCTION MANUAL

MANUEL D'INSTRUCTIONS

MANUAL DE INSTRUCCIONES

CARATTERISTICHE TECNICHE/TECHNISCHE MERKMALE/TECHNICAL CHARACTERISTICS/CARACTERISTIQUES TECHNIQUES/CARACTERÍSTICAS TÉCNICAS

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Piston-cylinder electrical compressor lubricated for life.
Metal casing painted with epoxy powder paint.

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The MINIASPEED BATTERY PRO aspirator is a professional portable aspirator, specific for secretion removal. It has a vacuum regulator (11), a vacuum gauge (12) and 1000 ml container (8) with a protection device to prevent liquid from entering the suction pump and interrupting the aspiration flow. Lifetime lubrication, easy to handle, simple to use, reliable, strong and silent. This model has a rechargeable battery and can also operate using a cigar lighter socket directly in a car or a boat.

The MINIASPEED BATTERY PRO aspirator has the following 3A accessories: a 1000 ml container with protection device, battery charger, cigar lighter power supply cord, multi-voltage power supply unit (miniaspeed battery plus), 130 cm long silicon connection tube that can be sterilised, **sterile, disposable** cannula, **sterile, disposable** manual flow regulator, **disposable** antibacterial filter, **disposable** bag.
N.B.: Only use original 3A accessories.

IMPORTANT WARNINGS

This is a medical device for use on medical prescription. It must be operated as indicated in this user instruction manual. It is important for patients to read and understand the information concerning using and maintaining the unit. If you have any questions, contact your stockist. MICROBIAL CONTAMINATION: in the case of illnesses with a risk of infection or microbial contamination, we advise each individual to use his own accessories. (Consult your doctor).

The manufacturer has made every effort to ensure that all its products are of the highest quality and are reliable. Nevertheless, as for all electrical appliances, fundamental safety standards must be observed:

- Children and persons who are not self-sufficient may only use the unit under the strict supervision of a responsible adult who has read this manual.
- Never use adapters for supply voltages different to the voltage shown on the data plate on the back of the unit. Keep the cord away from hot surfaces.
- The unit should not be used in the presence of inflammable anaesthetic mixtures with air, oxygen or nitrous oxide.
-  Never handle the supply cord plug with wet hands or use the unit when taking a bath or a shower. Never leave the unit close to water, do not submerge it in any liquids. If by chance it falls into water, pull the plug out of the socket immediately before recovering it. Do not use the unit if the plug or supply cord are worn or wet (send it immediately to your stockist).
- Unit casing not protected against liquid penetration.
- Only authorised personnel may perform maintenance and/or repair work. Unauthorised repairs annul the guarantee.
- Ensure that the connections and container closing are performed carefully to prevent suction losses.
- Do not tip the container over while connected to the unit when operating, as liquid may be sucked inside the appliance causing pump damage. If this happens, switch the aspirator off immediately and empty and clean the container. Send the appliance to your stockist.
- When the protection device intervenes aspiration is interrupted; empty the container (8) and perform the cleaning operations.
- The cannula and the manual control of the aspirated flow are **sterile, disposable** products: They must be replaced after every application.
- Check the use-by-date on the original packaging of the cannula and check the integrity of the sterile packaging.
- The **disposable** antibacterial filter must be replaced after every application.
- The power supply cord with cigar lighter plug has a safety fuse for inspection in the event of a fault.
- Never use the battery charger with other appliances or for uses other than that established by this manual. Never use MINIASPEED BATTERY PRO with other power supply units.
- ⚠ the appliance contain a hermetic lead battery. When disposing of the battery, follow the regulations in force for disposing of toxic-noxious refuse.

USE INSTRUCTIONS

Each time before use, ensure that all the accessories are perfectly clean according to the instructions indicated in the "CLEANING AND DISINFECTION OPERATIONS"

1. Operation with the cord and cigar lighter plug:

- 1.1 For use in a car, camper, boat, etc....., connect the power supply cord (1) to the appliance socket (15).
- 1.2 Before using the appliance, check the charge status of the vehicle battery.
- 1.3 Connect the appliance as indicated in figure 1.
- 1.4 By the flow regulator (11) turning the knob to "+" we obtain a higher vacuum and by turning it to "-" a lower one (bar/KPa); these values are read on the vacuum gauge (12).

- 1.5 Start the appliance by placing the switch to the "I" (ON) position (10) (green pilot light on).
1.6 After the application, switch off the appliance, detach the supply cord from the current socket (1) and perform the cleaning operations as illustrated in the "CLEANING AND DISINFECTION" paragraph
1.7 In the event of a blown fuse in the cigar lighter plug, replace with a 6A – 250V delayed fuse – dimensions Ø 6.3 x 30 mm, by unscrewing the end of the cigar lighter plug.

2. Operation with a battery and/or with the multi-voltage power supply unit:

2.1 The appliance is supplied with the battery partially charged. We therefore recommend charging it before use.

2.2 The battery is charged by connecting the multi-voltage power supply unit (3) to the mains supply and to the appliance socket (15), with the switch (10) switched to the (O) (OFF) position. Maximum charging time 6 hours with around 45 minutes of autonomy.

2.3 Operation only with the internal battery: if the green pilot light goes out during use, and the yellow pilot light comes on, this means that the remaining autonomy is around 10 minutes (reserve). Therefore terminate the application if possible, paying great attention to the flashing of the red pilot light. In this case do not use the appliance for more than 3 minutes to avoid damaging the battery. If the battery is flat and you want to continue the application, if possible, connect the multi-voltage power supply unit (3) to the appliance socket (15) and place the switch (10) to the "I" position.

2.4 At the end of the application, place the switch (10) to the "O" (OFF) position and leave the power supply unit (3) connected to the appliance socket (15) to charge the battery. When the appliance is not in use, we recommend leaving the power supply unit (3) connected so that the battery charge level is always optimum.

2.5 For aspiration operations see points 1.3; 1.4; 1.5.

The aspirator switches off automatically when the voltage falls below 9.5V to protect battery efficiency. In the event of extreme need, the operator can re-start the aspirator by pressing the start button which will remain on for 1 min.

4. Secretion collection vessel - 1000cc. (Fig. 2)

The 1000 cc. collection bottle supplied with the aspirator can be used in two ways: as a collection vessel which can be sterilised as shown in photograph 1A or as a collection vessel with disposable bag as illustrated in photograph 1B.

Sterilisable secretion collection vessel

The antibacterial filter (6) should be inserted directly into the vessel top. Do not use the aspirator without the antibacterial filter, because from a bacteriological point of view, it becomes dangerous for the patient. The collection vessel is supplied with an overflow valve, polypropylene plug (2) and polycarbonate vessel (1). The plug (2) foresees the direct insertion of the antibacterial filter, which can only take place on the connector called "VACUUM". The antibacterial filter also performs the function of protecting the aspiration circuit from any contaminating agents sucked in during use. All the components of the vessel can be sterilised using a conventional system in an autoclave at a temperature of 121°C, or by boiling for 10 minutes. We recommend replacing the complete vessel at every 30 sterilisation cycles. Do not overturn the vessel during use, in order to prevent the intervention of the non-return valve (3-4-5); should this occur, switch the aspirator off and detach the tube connected to the antibacterial filter. Never use the aspirator without the secretion collection vessel and/or without the antibacterial filter.

Connection: connect one end of the short, sterilisable silicon tube (10) to the antibacterial filter connector (6) and insert this latter into the "VACUUM" hole of the blue top (2); connect the other end to the "INLET" connector of the aspirator. Connect one end of the long sterilisable silicon tube to the "PATIENT" connector of the blue top (2); to the other end, connect the disposable, sterile manual regulator (9) and connect the disposable, sterile cannula (8) to this latter.

Secretion collection vessel with disposable bag. The aspirator can be used with the re-usable 1000 cc. transparent secretion collection vessel (1) and with the disposable bag supplied. In this case, the antibacterial filter is integrated in the disposable bag and therefore the antibacterial filter (6) and blue top with valve (2) are not used. The filter integrated in the bag also performs the function of preventing the return of aspirated liquids towards the aspirator. In this case, to restore device operation, it is necessary to replace the disposable bag. For the cleaning and disinfecting operations of the tubes and vessel, sterilise the single parts in an autoclave at a maximum temperature of 121°C, or by boiling for 10 minutes. The bag is disposable and must always be replaced after every use. The bag must be completely inserted in the vessel in order to prevent any vacuum losses.

N.B.: maximum disposable bag usage vacuum: -0.75 bar.

Connection: connect one end of the short sterilisable silicon tube (10) to the yellow connector (VACUUM) of the top (photograph 1B) and the other end to the "INLET" connector of the aspirator. Connect one end of the long sterilisable silicon tube to the red connector (PATIENT) and connect the disposable sterile manual regulator (9) and the disposable sterile cannula (8) to the other. **N.B.: only use Meditea 1 I. disposable bags - code**

M043002/A.

CLEANING AND DISINFECTION OPERATIONS

N.B.: If using chemical disinfectants, follow the manufacturer's instructions exactly.

- The cannula and the aspirated flow manual command are sterile, disposable products and must be replaced after every application.
- The antibacterial disposable filter must be replaced after every application.
- Never leave the appliance in water or submerged; clean the external casing of the appliance using only a damp cloth with detergent (non abrasive).

Electromagnetic compatibility – additional notes in conformity with the IEC 60601-1-2/A1:2004-09 Standard.

Emission		
Emission test	Conformità	Electromagnetic Environment - guidance
RF Emissions Cispr 11	Group 1	MINIASPEED BATTERY PRO uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions Cispr 14	--	MINIASPEED BATTERY PRO is suitable for use in all establishments and those directly connected to the public low-voltage power supply network that supplies buildings used to domestic purposes. It is possible to use the device in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, adopting provisions in the installation, such as a longer distance from potentially sensitive appliances.
Harmonic Emissions IEC 61000-3-2	Classe A Conforms	It is possible to use the device in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings user for domestic purposes.
Emissioni di fluttuazioni di tensione/flicker IEC 61000-3-3	Conforme	

Immunity aspects

MINIASPEED BATTERY PRO is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity test	Test level EN 60601-1-2	Compliance Level	Electromagnetic environment - guide
Electrostatic discharge (ESD) EN 61000-4-2	± 6kV contact ± 8kV air	± 6kV contact ± 8kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Burst/Fast Transient EN 61000-4-4	±2kV power supply lines	±2kV power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN 61000-4-5	±1kV differential mode ±2kV common mode	±1kV differential mode ±2kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines EN 61000-4-11	< 5% UT (>95% dip in UT) for 0,5 cycles 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (>95% dip in UT) for 5 seconds	< 5% UT (>95% dip in UT) for 0,5 cycles 40% UT (60% dip in UT) for 25 cycles 70% UT (30% dip in UT) for 25 cycles < 5% UT (>95% dip in UT) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions. It is recommended that the Laser/Pulsed Light BiosYag be powered from an uninterruptible power supply or a battery.
Power frequency magnetic field EN 61000-4-8	3 A/m	3 A/m	Magnetic power frequency fields should be that of a typical commercial or hospital environment.

Immunity aspects at r.f.

MINIASPEED BATTERY PRO is intended for use in the electromagnetic environment specified below. The customer or the user of the navigator should assure that it is used in such an electromagnetic environment.

Immunity test	Test Level EN 60601-1-2	Compliance Level	Electromagnetic environment - guide
RF conducted EN 61000-4-6	3 Veff from 150kHz to 80MHz	3 Veff from 150kHz to 80MHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from that equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1,2 \sqrt{P}$ da 150kHz a 80MHz $d = 1,2 \sqrt{P}$ da 80 MHz a 800 MHz $d = 2,3 \sqrt{P}$ da 800 MHz a 2,5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).
RF radiated EN 61000-4-3	3 V/m from 80MHz to 2,5GHz	3 V/m from 80MHz to 2,5GHz	

Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:

Recommended separation distances between portable and mobile RF communications equipment and the Laser/Pulsed Light Bios Yag			
<p>MINIASPEED BATTERY PRO is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.</p>			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of the transmitter (m)		
	From 150kHz to 80MHz $d = 1,2 \cdot \sqrt{P}$	From 80MHz to 800MHz $d = 1,2 \cdot \sqrt{P}$	From 800MHz to 2,5GHz $d = 2,3 \cdot \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer

Notes:
(1) At 80 MHz and 800 MHz the separation distance for the higher frequency range applies.
(2) These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

PROBLEMS, CAUSES AND SOLUTIONS

PROBLEMS	POSSIBLE CAUSES	SOLUTIONS
Excessive noise.	Damaged pump or blockages in the internal aspiration circuit	Send to the assistance circuit.
The unit switches on but does not aspirate	- Damaged pump - Vacuum regulator fully open. Connection tubes disconnected and/or badly connected, broken connection tubes. Container not in a vertical position, full, or defective overflow valve. Possible blockage of the hydraulic circuit inside the unit.	- Send to the assistance circuit. - Check the position of the vacuum regulator. Check the connections and the integrity of the tubes. Position the container in a vertical position, check the overflow valve (blocked) and/or replace the silicon tubes.
The vacuum rate cannot be regulated	Damage to the internal hydraulic system or blockage of the connection tubes to the aspiration unit.	Send to the assistance circuit.
When the appliance is switched on, the protection fuse always trips	Pump probably damaged or in short-circuit.	Send to the assistance circuit.
The vacuum gauge does not work	Liquid penetrating the pneumatic circuit.	Send to the assistance circuit.

Note: if you experience faults or malfunctioning problems different to those listed above, always and exclusively contact authorised assistance centres.

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GB**MINIASPEED BATTERY PRO:**

- 1) Cord with cigar lighter plug
- 2) Disposable manual regulator
- 3) Multi-voltage power supply unit
- 4) Power supply cord
- 5) Cannula
- 6) Bag with shoulder strap
- 7) Protection device
- 8) 1000 ml container
- 9) Container plug
- 10) ON – OFF switch
- 11) Vacuum regulator
- 12) Vacuum gauge
- 13) Antibacterial filter
- 14) 24 cm silicon tube
- 15) Cigar lighter socket and battery charger

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