ES - EN

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## ENGLISH

### WARNING

• Never under any circumstances attempt to open or inspect the internal components or accessories of the machine. If an inspection becomes necessary, please contact your supplier or a suitably competent technician.

• Never use the unit in humid or poorly ventilated environments, such as wet rooms, bathrooms etc. to prevent damage to the electrical components of your unit.

• Never use the unit outdoors; exposure to the elements can cause damage.

• If the power cord was destroyed, please ask a professional to service and replace it to avoid hazards, or contact the dealer for service.

### RECYCLING

 Whenever any of our electrical products or their presentation packaging bear this symbol, it means that these products should not be treated as conventional household waste in Europe. To ensure proper treatment of this waste, dispose of it in accordance with local laws or as required for the disposal of electrical equipment. This initiative will help preserve save natural resources and improve standards of environmental protection in regard to the treatment of electrical waste.

### WELCOME!

Thank you very much for purchase. Please read the instructions carefully before operating the instrument. You can also ask a professional for guidance to make sure you are using the instrument correctly. We will accept no responsibility for any accidents that occur due to incorrect operation.

The company reserves the right to amend the details of the instruments with no prior notice. If any mistakes are found in these instructions, please feel free to send us an e-mail informing us so that we can correct them.

### PACKAGE CONTENTS

• These components are subject to change without prior notice.

• Appearance of product may differ from the image.

### INSTRUCTIONS OF USE

**1.** Uncover the transparent lid and place the hair removal wax into the inner housing.

**2.** Plug the device to the main supply.

**3.** Turn the controller clockwise and the light will turn on.

**4.** The machine will heat the wax for about 20-30 min. When the wax reaches the melting point, the light will turn off. Then set the button to MID position.

**5.** Caution: Before starting the treatment, make sure that the wax temperature is safe. Let it cool down a bit or place the wax on the hair-removal paper band for treating.

**6.** Ensure that the skin is absolutely clean before its use.

**7.** Use a wooden spatula to cover the skin surface. The thinner the wax layer is, the easier the hair removal becomes.

**8.** Place the paper band onto the skin surface and pull it off quickly in the opposite direction of the hair growth.

**9.** After treatment, use a smoothing lotion onto the skin and massage.

### CAUTIONS

• Do not use on sensitive skins.

• Do not turn on the wax machine, if the wax pot is empty.

• When you add wax, move or clean the machine, turn it off.

• Switch off the power before unplugging the machine.

• The machine can only be cleaned with alcohols or specialised wax removal products.

• Do not leave children unattended with this appliance

© 2018 Weelko

Visit [www.weelko.com](http://www.weelko.com) for more information

Importer : Weelko Barcelona, S.L.

ESB65397556

### TECHNICAL SHEET

Weights and sizes	Packing size (cm)	17x17x17
	Gross weight (kg)	0,75
	Product size (cm)	15x15x12,5
Working characteristics	Capacity	400 ml
	Recommended operating time	20 - 30 min
Electrical characteristics	Rated voltage	220V ~ 240V
	Frecuency	50Hz/60Hz
	Rated power	130 W

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