

MIMLED 1000

Compact LED Surgery lamp. MIMLED 1000 offers 100.000 lux at 1 meter.

HIGH COLOUR RENDERING

The high colour rendering ($R_a > 94\%$), describes the exceptional quality of white light, allows tissues to be distinguished with optimal precision.

LIGHT INTENSITY DIMMING

Adjustable light intensity integrated in the handgrip and on the head which can be regulated using the panel control.

ENERGY SAVING LIGHT

Low power consumption.
33W MIMLED 1000

AVAILABLE IN 3 OPTIONS

- Trolley Stand.
- Wall version.
- Ceiling version (up to 3.10 meter height).



MIMLED 1000

Technical Details


Lux at 1 meter	100.000 lux
Lamp lifetime	50.000 hours
Colour temperature	4.500 °K
Colour rendering	≥ 94%
Light field diameter	195 mm
Power supply	100/230 V - 50/60 Hz
Consumption	33 W



*Adjustable light intensity
integrated in the grip*



Ceiling version
Cod. ML1000C



Wall version
Cod. ML1000W



Trolley Stand version
Cod. ML1000FL

