

MIMLED 600

Compact LED Surgery lamp. MIMLED 600 offers 60.000 lux at 1 meter.

HIGH COLOUR RENDERING

The high colour rendering ($R_a > 94\%$), describes the exceptional quality of white light, allows tissues to be distinguished with optimal precision.

LIGHT INTENSITY DIMMING

Adjustable light intensity integrated in the handgrip and on the head which can be regulated using the panel control.

ENERGY SAVING LIGHT

Low power consumption.

25W MIMLED 600

AVAILABLE IN 3 OPTIONS

- Trolley Stand.
- Wall version.
- Ceiling version (up to 3.10 meter height).



MIMLED 600

Technical Details	
Lux at 1 meter	60.000 lux
Lamp lifetime	50.000 hours
Colour temperature	4.500 °K
Colour rendering	≥ 94%
Light field diameter	195 mm
Power supply	100/230 V - 50/60 Hz
Consumption	25 W



Adjustable light intensity
integrated in the grip



Ceiling version
Cod. ML600C



Wall version
Cod. ML600W



Trolley Stand version
Cod. ML600FL