



LG EVA-CHAIR

Chiudibile / Closable

LG EVA-CHAIR è utilizzato principalmente per trasportare feriti o persone inabili in luoghi particolarmente stretti.

Realizzato in lega di alluminio ad alta resistenza, leggero, di piccole dimensioni;

Può essere piegato a piccolo, facilmente immagazzinabile;

4 ruote che facilitano il trasporto del paziente sul piano;

Con due manici pieghevoli sul retro e due sul davanti. La lunghezza delle maniglie anteriori in basso può essere regolata. /

LG EVA-CHAIR is mainly used to transport wounded or impaired people in particularly tight places.

Built with high resistance, lightweight aluminium alloy, compact size;

It can be folded up small, for easy storage;

4 wheels that make it easier to transport the patient on flat ground;

With two folding handles on the back and two at the front. The length of the bottom front handles is adjustable

LG EVA-CHAIR 37

NEW

La versione ridotta nei suoi 37 cm di larghezza è richiesta negli aeroporti come sedia da transito per la gestione dei pazienti in aereo. / The smaller version with a width of 37 cm is required for airports as a transit wheelchair to manage patients in aircrafts



LG EVA-CHAIR - Dati Tecnici / Technical Specifications

Ingombro max - Max overall dimensions	70x52x85 cm
Ingombro minimo (chiusa) - Minimum overall dimensions (closed)	90x54x20 cm
Peso netto - Net weight	11 Kg
Peso lordo - Gross weight	12 kg
Peso trasportabile - Transportable weight	160 kg

LG EVA-CHAIR 37 - Dati Tecnici / Technical Specifications

Ingombro "maniglie posteriori chiuse" / Size unfolded (LxDxH) "Back handles closed"	(74÷94,5) x 37 x 93,5 cm
Ingombro "maniglie posteriori aperte" / Size unfolded (LxDxH) "Back handles open"	(82,5÷103) x 37 x 93,5 cm
Ingombro minimo (chiusa) / Size folded (LxDxH)	19 x 37 x 93,5 cm
Peso netto / Net weight	10,4 kg
Peso lordo / Gross weight	12 kg
Peso trasportabile / Max load	160 kg