DRUGS & MEDICINE TROLLEYS - MADE IN ITALY







- 27494 PHARMACY TROLLEY small with 15 3-compartment partition
- 27495 PHARMACY-TROLLEY standard

with 20 3-compartment partition

• 27497 PHARMACY TROLLEY - large

with 30 3-compartment partition

A portable pharmacy, allowing to manage daily or weekly treatment of patients, with 15/20/30 3 compartments small drawers for easy drug storage (45-60-90 spaces).

Painted steel frame, technopolimer surface and panels. Can be completely closed by safety key lock.

Provided with 4 revolving castors Ø 125 (2 with brake); Size 820/1150x640xh 1,000 mm.

• 27498 PARTITION -

3 compartments - spare




DRUGS & MEDICINE CABINETS - MADE IN ITALY





- made in bilaminated board (27944/5) or melamine board (27942) 20 mm thick
- epoxy powder painted steel base frame with 4 adjustable feet Ø 50 mm
- 2 external lockable doors turning through 180° with 8 medicine compartments each (extendable to 10 compartments)
- perimetrical anti-dust gaskets
- inside 2 containers with 8 medicine compartments each side (extendable to 10 compartments) turning through 90° - dangerous drugs compartment provided with lock

GIMA code	MEDICINE CABINET					Material					Colour	
27942	Medicine cabinet					melamine particle board					white	
27944	Medicine cabinet					bi-laminated board					white	
27945	Medicine cabinet					bi-laminated board					any colour*	
27946	Compartment for medicine											
Colours:	412	431	441	452	2	482	484	499	854	861	1818	









- 27941 COMBO MEDICINE CABINETS With 24 medicine compartments one side and 6 ISO tray 1 shelf and 3 drawers
- made in bilaminated board 20 mm thick
- epoxy powder painted steel base frame with 4 adjustable feet Ø 50 mm Right side
- 1 external lockable door turning through 180° with 8 medicine compartments each (extendable to 10 compartments)
- perimetrical anti-dust gaskets
- inside 1 container with 8 medicine compartments each side (extendable to 10 compartments) turning through 90° provided with lock

Left side

- glass door, internal side with runners for ISO trays
- 3 drawers on the bottom Supplied with 6 ISO trays different heights (2xh 5 cm, 3xh 10 cm and 1xh 20 cm) and 1 shelf.
- 27946 COMPARTMENT FOR MEDICINE
- spare