24

FIT

PROFESSIONAL ALCOHOL TESTER

• 24491 FIT-333 ALCOHOL TESTER with printer

- 24486 AT500 EASY RECALIBRATION ALCOHOL TESTER with printer available jan 2018
- 24487 PRE-CALIBRATED SENSOR MODULE for 24486 spare
- 24494 BLOWING TUBES for 24486, 24491 box of 100 pcs.

Size: 188x78x33 mm Weight: 300 g

• 24495 PAPER ROLLS for 24486, 24491 - box of 10 Professional fuel cell sensor breath alcohol tester with printer and touch screen. The fuel cell is a porous disk with a thin layer of platinum. Any alcohol that comes in contact with the cell's surface is broken down chemically. For each broken down molecule of alcohol a given number of electrons are freed.

The numbers of electrons generated are in proportion to the amount of alcohol. Self recalibration (only 24486): easy replacement of pre-calibrated module.



24491 24486 Colour 3.2 TFT touch screen LCD Colour 3.5" touch screen LCD 5 keys operation: "Up", "Down", "Left", "Right", "OK" Navigation key input with touch panel Rechargeable AA battery (about 2500 mAh) Rechargeable 7.4 V Li-ion battery

Chamber & Sensor warm up automatically Warm up: Auto flow monitoring: Automatic flow detection Operating environment:

-20° to +70° (optimal -20° to +40°) -20 to +45°C (optimal 5 to 40°C) USB port -Data download: Windows XP, 7** -

10,000-18,000 memories for approx 3,000 full tests 50.000 memories full test

driver name, vehicle number, inspector name, test result with print date and time, location, calibration data, Measurement precision: Standard deviation ±0,025 mg/L or ±5% of the measurement

— Drift: typically <0,8% of the measurement value/month —</p> Calibration*: by appointed service by replacement of pre-calibrated sensor module

Unit setting: mg/L, mg/100 L, %BAC, %BAC *Need calibration every 12 months **Can export files as Microsoft Excel documents directly (24486 only)

Alcohol testers1 are supplied with integrated printer, 25 mouth pieces², 1 large mouth piece3, power supply4, power supply 100-240V, 50/60 Hz5, USB car adaptor⁷, wireless transmission module for FIT-333 printer^{6, 9}, carrying pouch¹⁰, 4 rolls of paper¹¹, charger¹².









• 24500 ALCOHOL TESTER

Display:

Keyboard:

Power Supply:

Compatibility:

Memories:

Compact device to measure your alcohol concentration in your breath.

Digital/graphic interface and animation, one-hand operation.

% BAC and mg/l values.

Supplied in carton box with 3 blowing tubes Ø 7 mm and one cleaning cap. Multilanguage box and manual: GB, FR, IT, ES, PT

- 24501 BLOWING TUBE box of 12 pcs.
- 32200 "AAA" BATTERIES blister of 4 pieces
- 24503 COUNTER DISPLAY Italian
- 24504 COUNTER DISPLAY English

SELF MONITORING

Large display to read exact value of your breath (3 colour backlight)



Counter display fits 6 tester

TECHNICAL SPECIFICATIONS

Range: Low & 0.15 ~ 1.50 % BAC or 0.075 ~ 0.75 mg/l Accuracy: ±0.2% BAC full range under alert level Battery: 3 x 1.5V AAA batteries (not included) Display: Graphic & Digital display with 3 decimal points (eg. 0.025%BAC)

Tolerance compensation: Pre-alert function together with tri-colour backlighting

Breath sampling: less than 5 seconds continuous breath sample

Recycle time: about 20 seconds Size: 107x45x25 mm approx. Weight: 50 g (without batteries) Alert level calibrated: 0.5 % BAC or 0.25mg/l



POCKET ALCOHOL TESTER

24499 POCKET ALCOHOL TESTER Pocket, easy to use alcohol tester

TECHNICAL SPECIFICATIONS

Range: Low & 0.2 ~ 1.5 %BAC Accuracy: ± 0.2% BAC at the calibrated point (alert level)

Power input: 2 x 1.5V AAA batteries (not included) Display: Digital display

Tolerance compensation: pre-alert function when close to alcohol limit

Auto power off: 2 minutes Breath sampling with auto flow detect: about 6 seconds continuous breath sample

Cycle time: about 20 seconds Alert level calibrated: 0.5% BAC

POCKET, ONLY 8.5x3x1 cm

SALIVA ALCOHOL STRIPS

24502 SALIVA ALCOHOL TEST STRIPS - box of 25 strips

GIMA Saliva Alcohol Test Strip is a rapid highly sensitive test for the determination of relative blood alcohol concentration in saliva specimens.

- results in just 3-5 minutes

- colour chart provides approximate relative blood alcohol concentrations (BAC) from 0.02% to 0.30%

- no instrumentation is required
- saliva testing avoids urine collection
- short expiry date (average 10-22 months)
- professional use





DRUG OF ABUSE SALIVA MULTITEST - 8 DRUGS

 ϵ

• 24566 MULTI DRUG SALIVA MIDSTREAM TEST oral fluid (saliva) - box of 20 tests

A multi parameters, all-in-one test for simultaneous detection

DRUG OF ABUSE SALIVA MULTITEST Drug Cut off Cut off Drug Opiates 40 ng/ml 30 ng/ml Cocaine Cannabinoids (THC) Morphine 40 na/ml 25 ng/ml 6-MAM 4 na/ml Amphetamine, m-Amphetamine 40 ng/ml Hydrocodone and Oxycodone 40 ng/ml MDMA, MDA 50 ng/ml

of 8 different drugs in 10 minutes. Cost effective system for the screening of workers for illicit substance consumption.

