

# PROFESSIONAL ASPIRATORS - HIGH VACUUM



28280-1

CE 0434  
EN ISO 10079-1  
Class I  
Type B



28244-6

1 l Autoclavable Jar

## ASPEED PROFESSIONAL SURGICAL ASPIRATORS

Piston-type continuous cycle electric aspirators give high performance and great durability. Equipped with a protective thermal cut-out relay. They require no maintenance or lubrication. A motor-protector cap totally prevents aspirated

liquids or secretions from reaching and damaging the vacuum pump. Ideal for clinical use. Made in Italy.

GIMA code	ASPEED ASPIRATORS	Power	Pump	Case
28244	Aspeed 15 l	230 V	single	metal
28245	Aspeed 22 l	230 V	double	metal
28246	Aspeed 22 l	110 V	double	metal
28280	Aspeed 2 15 l	230 V	single	plastic
28281	Aspeed 2 22 l	230 V	double	plastic

GIMA code	ACCESSORIES AND SPARES FOR ASPEED
28255	Hydrophobic 99% bacterial filter
28258	Autoclavable jar 1l with cover
28259	Disposable suction liner
25480	Silicone tube 6x12 mm - roll of 30 m
28254	Adaptor to connect catheter to silicone tube
28261	Adaptor kit to use jar 1l (28258) with ASPEED manufactured before 2007

**STANDARD ACCESSORIES**  
Autoclavable polycarbonate jar 1,000 cc with safety valve (overflow protection)  
Disposable suction liner 1 l  
99% Antibacterial hydrophobic filter  
Sterile disposable cannula  
Sterile manual flow regulator  
Set of atoxic sterilizable silicone tubes  
Power Cable  
User Manual GB, FR, IT, DE, ES

	ASPEED		ASPEED 2	
	28244	28245/6*	28280	28281
Operating voltage:	230 V-50 Hz		230 V - 50/60 Hz	
	110-60 Hz		other voltage on request	
Bottle capacity:	1 l	1 l	1 l	1 l
High vacuum:	low flow	low flow	low flow	high flow
Air flow:	15 l/min	22 l/min	15 l/min	22 l/min
Adjustable vacuum level:	0 ÷ -0.85 bar (0 ÷ -85 kPa)			
Weight:	3.5 kg	4.5 kg	2.5 kg	3.2 kg
Case material:	metal		plastic	
Noise level:	55 dBA	65 dBA	55 dBA	55 dBA



# SUCTION ASPIRATORS HIGH VACUUM, LOW AND HIGH FLOW

- 28222 TOBI - suction aspirator  
220-230 V - 50/60 Hz
- 28224 SUPER TOBI - suction aspirator  
220-230 V - 50/60 Hz

Portable suction aspirators, ideal for tracheotomy and small surgery. Vacuum continuously adjustable with vacuum indicator. Have unbreakable 1,000 ml standard bottle (2,000 ml optional) autoclavable at 120°C with safety float control valve to prevent overflow. Silicone connection tube. ABS plastic case. Made in Italy



Silicone tube 6x10 mm (see page 159)  
Adaptor for catheter (code 28252)

STANDARD ACCESSORIES	Tobi		Super Tobi	
	Tobi	Super Tobi	Tobi	Super Tobi
Bottle 1,000 ml with cover	1	1	1	1
Antibacterial Filter	1	1	1	1
Suction catheter	1	1	1	1
Silicon Tube set	1	1	1	1
User manual GB, FR, IT, DE, ES	1	1	1	1

TECHNICAL SPECIFICATIONS	28222 Tobi		28224 Super Tobi	
	Operating voltage:	220-230 V - 50/60 Hz		220-230 V - 50/60 Hz
Power consumption:	184 W		106 W	
Bottle capacity:	1,000 ml		1,000 ml	
Flow (air litres/min):	18 l/min		40 l/min	
Max suction:	-0.75 bar (563mm/Hg)		-0.80 bar (600mm/Hg)	
Working time: minutes	20 ON / 40 OFF		120 ON / 60 OFF	
Size (cm):	37x22xh 21 cm		37x22xh 21 cm	
Weight:	3.5 kg		5.25 kg	
Norms:	IEC 601-1		IEC 601-1	