HOT AIR STERILIZERS: GIMETTE 1.5, 7, 21, 28 AND 50 LITRES



• 35601 GIMETTE 1.5 - hot air sterilizer • 35607 GIMETTE 7 - hot air sterilizer

Automatic Gimette are compact and high quality hot air sterilizers with the same reliability and efficiency as larger models. Outside structure is made of painted plate with anti-scratch epoxy powder. Stainless steel lid and sterilization chamber with approx. Capacity of 1.5 I and 7 I.

Removable, perforated for air circulation, tray with pliers. All components certified by most important European brands, assembling in compliance with EN-61010. A second safety thermostat triggers in case temperature should exceed 220°C. Red LED indicates any anomaly in thermoregulation and consequent starting of safety thermostat.

- 35631 RESISTOR 250 W for 35601 spare
- 35632 RESISTOR 400 W for 35607 spare
- 35612 GASKET for 35607 spare
- 35613 GASKET for 35601 spare



TECHNICAL SPECIFICATIONS

GIMETTE 1.5 245x120x50 mm = vol. 1.5 dm3 Inside dimensions:

Outside dimensions: 310x180x180 mm Weight: 4.3 kg Power consumption: 250 Watt Fuse: 1 A

GIMETTE 7 Inside dimensions: 410x180x90 mm = vol. 7 dm3

Outside dimensions: 490x310x180 mm Weight: 8 kg Fuses: 2 x 2.5 A Power consumption: 400 Watt

high precision with capillary bulb Same specifications for both models:

Operating Voltage: 230 V - 50 Hz

automatic, adjustable 0-200 °C Thermostat: Timer: manual from 0 to 120 minutes

2° Thermostat: bi-metallic certified

Green LED: timer, orange LED: heating, red LED: alarm

Made in Italy - Warranty 1 year User/Service manual in IT, GB, FR

Thermometer:



- 35621 GIMETTE 21 hot air sterilizer
- 35628 GIMETTE 28 hot air sterilizer
- 35629 GIMETTE 50 hot air sterilizer

Fully automatic hot air sterilizer, with adjustable temperature from 0°C to 200°C, kills off even particularly resistant micro-

C E 0425

Stainless steel structure is powder painted, except for the door, to assure maximum heat and shock resistance over time. Four perforated, anodized aluminium trays make hot air circulation easier. Capacity approximately 20 I, 28 I and 50 I. A second safety thermostat triggers in case temperature should exceed 230°C.

SPARE ACCESSORIES			
STERILIZER	TRAY	GASKET	RESISTOR 600 W
Gimette 21	35610	35614	35639
Gimette 28	35611	35615	35639
Gimette 50*	-	35617	35639

*Needs 2 resistors



Red LED indicates any anomaly in thermoregulation and consequent starting of safety thermostat.

TECHNICAL SPECIFICATIONS

GIMETTE 21 GIMETTE 28 GIMETTE 50 300x250x260 380x250x260 405x340x350 Inside dimensions mm: Approx. volume 20 dm3 28 dm³ 50 dm3 440x340x400 605x460x480 Outside dimensions mm: 550x340x400 Weight: 13 kg 17 kg 24.5 kg Power consumption: 600 Watt 500 Watt 1200 Watt

Same specifications for 3 models: 230 V - 50 Hz Operating Voltage:

2 x 4 A (35621, 35628) 2 x 8 A (35629) Fuses:

manual from 0 to 120 minutes Timer:

Thermostat: automatic

Thermometer: high precision with capillary bulb 2° Thermostat: bi-metallic certified

Green LED: timer, orange LED: heating, red LED: alarm

Made in Italy - Warranty 1 year User/Service manual in IT, GB, FR