

## SKIN MARKER



**• 33179 SKIN MARKER - dual tips**  
**gentian violet - box of 10**  
 Double tips skin marker with polypropylene cap and body. Polyester fiber tips.  
**Ultra fine tip utility marker** to write on adhesive labels and infusion bags.

**Regular tip surgical skin marker** for use on patient skin.  
 Washable with water and soap.  
 Size 138xØ 12.5 mm.  
 Multilanguage box and manual: GB, FR, IT, DE, ES, PT, GR, PL, SE, RO, HU, Arab.

## DERMATOGRAPH MARKERS




- 33175 DERMATOGRAPH MARKER - blue
- 33176 DERMATOGRAPH MARKER - white
- 33177 DERMATOGRAPH MARKER - red
- 33178 DERMATOGRAPH MARKER - mix\*

\*2 pencils each colour

Special dermatograph pencil to mark the skin for identification of areas to be treated. Made in cedar wood with

wax lead. Size 175xØ 7 mm. Suitable for general surgery, dermatology, orthopaedy, radiology, ...  
 Pencils are sold non sharpened in box of 6 with 1 pencil sharpener included.  
 Multilanguage box and manual: GB, FR, IT, DE, ES, PT, GR, PL, SE, RO, HU, Arab.

## SKIN MARKER - STERILE



**• 33181 CLINICON® - fine tip - sterile - box of 10**  
 The Clinicon® Surgical Skin Marker features

a versatile tip for precise fine line marking preferred in plastic surgery or when held at an angle, broad line marking for applications in general, neuro and vascular surgery, as well as radiology and radiation therapy. The marker's high visibility Gentian violet ink is nontoxic, non tattooing and surgical prep-resistant. Markings can be removed from the skin with isopropyl (rubbing) alcohol. The package is easy to open, HIBC bar coded and includes a separate six inch, flexible ruler calibrated in both inches and centimeters. Sterile EO.

## ETHICON SUTURES - ABSORBABLE AND NON ABSORBABLE



# ETHICON



### NON ABSORBABLE STERILE SURGICAL SUTURES

| GIMA code  | ETHICON code | Gauge | Metric | Colour | Circle   | Needle length | Needle     | Suture length | Minimum order |
|--|--------------|-------|--------|--------|----------|---------------|------------|---------------|---------------|
| <b>ETHILON - monofilament suture</b>                 |              |       |        |        |          |               |            |               |               |
| 22300  | 660H         | 6/0   | 0.7    | Black  | 3/8      | 16 mm         | FS-3       | 45 cm         | box of 36     |
| 22301  | 661H         | 5/0   | 1.0    | Black  | 3/8      | 19 mm         | FS-2       | 45 cm         | box of 36     |
| 22302  | 662H         | 4/0   | 1.5    | Black  | 3/8      | 19 mm         | FS-2S      | 45 cm         | box of 36     |
| 22303  | 663H         | 3/0   | 2.0    | Black  | 3/8      | 24 mm         | FS-1       | 45 cm         | box of 36     |
| 22304  | EH7506H      | 3/0   | 2.0    | Black  | 3/8      | 19 mm         | FS-2       | 45 cm         | box of 36     |
| 22305  | EW661        | 5/0   | 1.0    | Blue   | 3/8      | 19 mm         | FS-2       | 75 cm         | box of 36     |
| <b>PERMA-HAND - silk sutures - braided</b>           |              |       |        |        |          |               |            |               |               |
| 22310  | 18501G       | 4/0   | 1.5    | Black  | 1/2      | 17 mm         | J-1        | 45 cm         | box of 12     |
| 22311  | 18502G       | 3/0   | 2.0    | Black  | 1/2      | 22 mm         | X-1        | 45 cm         | box of 12     |
| 22312  | 18505G       | 2/0   | 3.0    | Black  | 1/2      | 22 mm         | X-1        | 45 cm         | box of 12     |
| 22313  | 18507G       | 4/0   | 1.5    | Black  | Straight | 19 mm         | ST-4       | 45 cm         | box of 12     |
| 22314  | 18508G       | 4/0   | 1.5    | Black  | 3/8      | 13 mm         | P-3 prime  | 45 cm         | box of 12     |
| 22315  | 18509G       | 3/0   | 2.0    | Black  | 1/2      | 22 mm         | SH-1 Plus  | 45 cm         | box of 12     |
| 22316  | EH7149G      | 4/0   | 1.5    | Black  | 3/8      | 19 mm         | FS-2       | 45 cm         | box of 12     |
| <b>PROLENE - blue pigmented monofilament sutures</b> |              |       |        |        |          |               |            |               |               |
| 22320  | EH7150H      | 4/0   | 1.5    | Blue   | 3/8      | 19 mm         | FS-2       | 75 cm         | box of 36     |
| 22321  | EH7151H      | 4/0   | 1.5    | Blue   | 3/8      | 19 mm         | FS-2       | 45 cm         | box of 36     |
| 22322  | W8873T       | 6/0   | 0.7    | Blue   | 3/8      | 16 mm         | PC-3 Prime | 45 cm         | box of 24     |

Non absorbable surgical sutures are indicated for use in general soft tissue approximation and/or ligation included use in cardiovascular, ophthalmic and neurosurgical procedures

### ABSORBABLE STERILE SURGICAL SUTURES

| GIMA code   | ETHICON code | Gauge | Metric | Colour | Circle | Needle length | Needle    | Suture length | Minimum order |
|---|--------------|-------|--------|--------|--------|---------------|-----------|---------------|---------------|
| <b>Monocryl is indicated for use in general soft tissue approximation and/or ligation, but not for cardio, neuro, micro and ophthalmic surgery</b>                                |              |       |        |        |        |               |           |               |               |
| <b>MONOCRYL - complete absorption in 91-119 days</b>  |              |       |        |        |        |               |           |               |               |
| 22350   | Y493H        | 5/0   | 1.0    | Undyed | 3/8    | 19 mm         | P-3 prime | 45 cm         | box of 36     |
| <b>Vicryl Plus is indicated for use in general tissue approximation and/or ligation except for ophthalmic, neuro, cardiovascular</b>  |              |       |        |        |        |               |           |               |               |
| <b>VICRYL PLUS BRAIDED - complete absorption in 56-70 days</b>  |              |       |        |        |        |               |           |               |               |
| 22360   | VCP392ZH     | 4/0   | 1.5    | Violet | 3/8    | 17 mm         | FS-2      | 45 cm         | box of 36     |
| 22361   | VCP305H      | 3/0   | 2.0    | Violet | 1/2    | 17 mm         | RB-1 Plus | 70 cm         | box of 36     |
| 22362   | VCP452H      | 2/0   | 2.0    | Violet | 3/8    | 24 mm         | FS-1      | 70 cm         | box of 36     |
| 22363   | VCP453H      | 2/0   | 3.0    | Violet | 3/8    | 24 mm         | FS-1      | 70 cm         | box of 36     |
| <b>Vicryl Rapid Braided indicated only for: use in superficial soft tissue approximation of the skin and mucosa, where only short term (7-10 days) wound support is required.</b> |              |       |        |        |        |               |           |               |               |
| <b>VICRYL RAPID BRAIDED - complete absorption in 42 days</b>  |              |       |        |        |        |               |           |               |               |
| 22370   | V76H         | 3/0   | 2.0    | Undyed | J-1    | 45 cm         | FS-2      | 45 cm         | box of 36     |
| 22371   | V2920G       | 4/0   | 1.5    | Undyed | 3/8    | 19 mm         | FS-2      | 45 cm         | box of 12     |
| 22372   | V2930G       | 3/0   | 2.0    | Undyed | 3/8    | 19 mm         | FS-2      | 45 cm         | box of 12     |
| 22373   | V2150G       | 3/0   | 2.0    | Undyed | 1/2    | 17 mm         | RB-1 Plus | 70 cm         | box of 12     |
| 22374   | V2190G       | 3/0   | 2.0    | Undyed | 1/2    | 22 mm         | SH-1 Plus | 70 cm         | box of 12     |
| 22375   | VR229H       | 3/0   | 2.0    | Undyed | 3/8    | 14 mm         | FS-2      | 75 cm         | box of 36     |