

mindray

Model NO. 0146-00-0091-01

Device Identifier:LI23I002A

DC 11.1V \approx 4800mAh 53.28Wh

可充电锂离子电池 Rechargeable Li-ion Battery

3ICR19/66-2



中国制造
MADE IN CHINA

警告：

使用不当可能导致电池起火、爆炸或燃烧。切勿挤压、刺破、拆卸或使电池短路。切勿放置于火中或水中。在-20°C (-4°F) 至45°C (113°F) 温度下存储。只能使用指定的充电器。请阅读本产品的操作说明。

CAUTION:

Risk of fire, explosion, or burns. Do not crush, puncture, disassemble or short circuit the battery. Do not dispose of the battery in fire or water. Store the battery in the -20°C (-4°F) to 45°C (113°F) environment. Use the specified charger.



深圳迈瑞生物医疗电子股份有限公司
SHENZHEN MINDRAY BIO-MEDICAL ELECTRONICS CO., LTD.

047-010822-02(1.0)