



## Innovation Deluxe Daycare



The Innovation Deluxe Daycare couch is designed for specialist medical procedures, critical care and oncology units. The fully electric functionality provides smooth lifting and patient profiling during patient examination and treatment. The couch features 2-way tilt and auto-electric CPR facility with auto-levelling software.

### Models

- MG3693** Electric profiling of height, backrest and foot section with auto-levelling, 2-way tilt and CPR facility. Large central locking wheel design with steering facility
- MG3493** Electric profiling of height, backrest and foot section with auto-levelling, 2-way tilt and CPR facility. Large independently braked wheel design



Model MG3693  
with large central locking wheel  
design with steering facility

### Deluxe Features

- Electric height range from 47cm to 97cm
- Electric trendelenburg tilt to  $\pm 15^\circ$
- Electric profiling backrest angle from  $0^\circ$  to  $+80^\circ$
- Electric footrest angle from  $+10^\circ$  to  $-35^\circ$
- Safe working load and lifting capacity of 250Kg (550lbs)
- Auto-electric CPR function via handset
- Auto-intelligent operation
- Vertical lifting design to minimise space and provide maximum stability throughout all profiling positions
- Smooth moulded base cover minimises dirt traps and facilitates easy cleaning
- Anti-bacterial powder coated frame and base cover
- Large Tente®  $\varnothing 125\text{mm}$  castors
- **SEERS ComfortPlus™** Upholstery - features anti-microbial vinyl incorporating pressure relieving design for additional patient comfort with a foam depth of 60mm
- Standard width upholstery of 70cm
- Clearance beneath frame for mobile hoist
- Head support cushion fitted on Deluxe models
- Paper roll holder fitted on Deluxe models
- Fold down side support arm rests for patient comfort and security
- Choice of 15 standard upholstery colours

### Optional Accessories

- 6054** Telescopic IV pole
- 6057** Battery back-up facility
- 6241** Patient lap tray
- 6242** Adjustable foot board

See page 48-49 for additional information on our range of optional accessories.



Model 6241 optional  
patient lap tray



Model 6242 optional  
adjustable foot board

2-way trendelenburg tilt  
facility to  $\pm 15^\circ$



Low patient transfer height