



Innovation Deluxe Dialysis



The Innovation Deluxe Dialysis couch has been designed for use within specialist renal dialysis environments. The couch features auto-electric CPR and fully electric profiling with a 2-way tilt facility to ensure a smooth transition when positioning patients for dialysis procedures and treatments.

Models

- MG3690** Electric profiling of height, backrest and foot section with auto-levelling, 2-way tilt and CPR facility. Large central locking wheel design with steering facility
- MG3490** Electric profiling of height, backrest and foot section with auto-levelling, 2-way tilt and CPR facility. Large independently braked wheel design



Model MG3490 with large independently braked wheel design

Deluxe Features

- Electric height range from 47cm to 97cm
- Electric trendelenburg tilt to $\pm 15^\circ$
- Electric profiling backrest angle from 0° to $+80^\circ$
- Electric footrest angle from $+10^\circ$ to -35°
- Safe working load and lifting capacity of 250Kg (550lbs)
- Auto-electric CPR function via handset
- Auto-intelligent operation
- Vertical lifting design to minimise space and provide maximum stability throughout all profiling positions
- Smooth moulded base cover minimises dirt traps and facilitates easy cleaning
- Anti-bacterial powder coated frame and base cover
- Large Tente® $\varnothing 125\text{mm}$ castors
- **SEERS ComfortPlus™** Upholstery - features anti-microbial vinyl incorporating pressure relieving design for additional patient comfort with a foam depth of 60mm
- Standard width upholstery of 70cm
- Clearance beneath frame for mobile hoist
- Head support cushion fitted on Deluxe models
- Paper roll holder fitted on Deluxe models
- Multi-adjustable arm rests for patient comfort
- Adjustable upholstered foot board for patient comfort
- Choice of upholstery colours, for full range of colour options see page 51

Optional Accessories

- 6054** Telescopic IV pole
- 6057** Battery back-up facility
- 6241** Patient lap tray

See page 48-49 for additional information on our range of optional accessories.



Low patient transfer height



Fully profiled for patient comfort



2-way trendelenburg tilt facility to $\pm 15^\circ$