



\*When fitted with 9400 Easy Clean Upholstery

## Innovation Standard 2 Section



The Innovation Standard 2 section couch features fully electric profiling of height and backrest angle. This couch has the ability to support very heavy patients and is suitable for general examinations and minor treatment procedures. The anti-bacterial base cover is supplied as standard on all Innovation models.

### Models

- MG2670** Electric height with electric backrest. Large central locking wheel design with steering facility
- MG2470** Electric height with electric backrest. Large independently braked wheel design



Model MG2470 with independent braked castors. Optional head support cushion and side support rails

### Standard Features

- Safe working load and lifting capacity of 250Kg (550lbs)
- Electric height range from 45cm to 95cm
- Electric profiling backrest angle from -20° to +85°
- Designed to provide maximum stability throughout all profiling positions
- Smooth moulded base cover minimises dirt traps and facilitates easy cleaning
- Anti-bacterial powder coated frame and base cover
- Anti-microbial, fire retardant upholstery
- Standard width upholstery of 65cm with 40mm depth medium density foam
- Large Tente® Ø125mm castors
- Option of independent braked wheel design or central locking wheel design with steering facility
- Clearance beneath frame for mobile hoist
- Choice of upholstery colours, for full range of colour options see page 51

### Optional Accessories

- 6003** Side support rails (Pair)
- 6003P** Side support rail cushions (Pair)
- 6004** Breathing hole and plug for prone treatments
- 6005** Paper roll holder
- 6006** Upper limb support
- 6007** Removable head/knee support
- 6054** Telescopic IV pole
- 6057** Battery back-up facility

See page 50-51 for additional information on our range of optional accessories



Low patient transfer height



Negative head section to -20°