



Medicare 3 Section Couches



The Medicare 3 Section Couch provides a very stable and versatile support surface for prone, supine and seated treatment procedures and examinations. Both the elevating sections can be set to a wide range of angles for specific patient positioning.

Models

- SM3550** Hydraulic with manual backrest
- SM3560** Electric with manual backrest
- SM3570** Electric with electric backrest



Model SM3560 Colour cream

Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Adjustable backrest from horizontal to +85°
- Adjustable angle head/foot section from -25° to +75°
- The power assisted backrest and foot sections incorporate anti-trap design on Models SM3550/SM3560
- Choice of manual or electric height elevation
- Excellent height range from 45cm to 98cm
- Retractable wheel design and adjustable foot for stability
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Frames pre-drilled for addition of accessories
- Anti-bacterial powder coated frame
- Standard width 65cm upholstery
- Choice of upholstery colours, for full range of colour options see page 51

Optional Accessories

- 6053** Protective base frame cover
- 6004** Breathing hole and plug for prone treatments
- 6001** Perimeter foot switch operation
- 6200** Single foot switch operation for Model SM3560
- 6207** Double foot switch operation for Model SM3570
- 6005** Paper roll holder
- 6003** Side support rails (Pair)
- 6003P** Side support rail cushions (Pair)
- 6006** Upper limb support
- 6007** Removable head/knee support
- 6057** Battery back-up facility
- 6054** Telescopic IV pole

See page 50-51 for additional information on our range of optional accessories



Negative head/ foot section to -25°



Low patient transfer height