

MC6158
45cm to 59cm

MC6159
54cm to 74cm



Backrest
25cm
38cm

Seat Pad
40cm
40cm

Surgeons & Sonographers Chair



The surgeons and sonographers chairs feature an ergonomically designed seat pad and smaller backrest that will not impede the movement of the arm during scans and procedures. The curvature design and seat pad adjustment encourages pelvic tilt, prevents slouching and promotes positive posture. The chair is standard with a small, compact base to reduce trip hazards and features a wipe clean upholstered underside.

Models

MC6158 Surgeons chair - standard model
MC6159 Surgeons chair - high model
Foot support ring fitted as standard



Fully upholstered backrest and seat pad for ease of cleaning



Model MC6158 Surgeons and Sonographers chair, standard height model



Model MC6159 Surgeons and Sonographers chair, high model with foot support ring



Model MC6159 Surgeons and Sonographers chair with optional foot height control and armrests

Features

- Wide choice of models available
- Ergonomically designed seat
- Compact backrest for improved scanning procedures
- Seat and backrest upholstered on both sides
- Weight capacity and safe working load of 150Kg
- Anti-bacterial and waterproof vinyl upholstery
- Restraint castors fitted as standard

Optional Accessories

- 6178** Anti-static castors
- 6455** Load locking castors
- 6172** Manual locking castors
- 6173** Glides for added stability on high models
- 6175** Height adjustable armrests
- 6174** Foot height control operation
- 6176** Foot support ring for standard model