



## Therapy 3 Section Drop End Couch



An extremely strong and rigid 3 section couch incorporating both an adjustable angle backrest and a single drop end head section, these models are suitable for a wide range of manual therapy and patient examination procedures. The couch has an excellent height range of 48cm to 101cm and is fitted with a breathing hole and plug as standard.

### Models

- ST3556** Hydraulic with manual backrest
- ST3566** Electric with manual backrest



Model ST3566 with electric elevation.  
Colour dark blue

### Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Adjustable angle backrest from horizontal to +90°
- Adjustable angle head section from +10° to -90°
- Power assisted backrest and head sections incorporate anti-trap design
- Choice of hydraulic or electric height elevation
- Excellent height range from 48cm to 101cm
- Retractable wheel design with adjustable foot for stability on uneven surfaces
- Breathing hole and plug included as standard
- Purpose made precision bearings on all pivot points
- Clearance beneath frame for mobile hoist
- Hand switch fitted as standard on electric models
- Standard width 65cm upholstery

### Optional Accessories

- 6001** Perimeter Foot Switch
- 6200** Single Foot Switch
- 6239** Dual Foot Switch Operation
- 6006** Upper Limb Support
- 6002** Side Support Rails
- 6002P** Side Support Rail Cushions
- 6005** Paper Roll Holder
- 6007** Removable Head/Knee Support
- 6053** Protective Base Shroud
- 6010** Couch Cover (Breathing Hole)
- 6011** Couch Cover (Without Breathing Hole)

See page 48-55 for additional information on our range of optional accessories and colour options



Negative head/foot section to -90°



Couch in horizontal position