



Therapy Bariatric 2 Section Couch



The Bariatric 2 Section Couch has a safe working load and lifting capacity of 325Kg (715lbs). Available with either a retractable wheel system or large central locking castors with steering facility. These models have a width of 80cm and the breathing hole and plug is fitted as standard.

Models

- ST2573** Electric with electric backrest
Retractable wheels
- ST2673** Electric with electric backrest
Large central locking castors with steering facility
- ST6207** Dual foot switch option



Model ST2573 Bariatric 2 section with retractable wheel system. Colour dark blue

Therapy Bariatric 3 Section Couch



The Bariatric 3 Section Couch features electric height operation, head/foot section and backrest angle adjustment with a safe working load and lifting capacity of 325Kg (715lbs). Available with either a retractable wheel system or large central locking castors with steering facility. These models have a width of 80cm and the breathing hole and plug comes fitted as standard.

Models

- ST3573** Electric with electric backrest and head/foot section. Retractable wheels
- ST3673** Electric with electric backrest and head/foot section. Large central locking castors with steering facility
- ST6209** Triple foot switch option



Central locking castors



Model ST3573 Bariatric 3 section with retractable wheel system. Colour light grey