



## Therapy Deluxe Couch



A very strong and rigid couch with full length top frame for added stability. These models are available with a choice of three head section designs with the option of split leg facility and have the perimeter foot switch operation supplied as standard.

### Models

<b>ST3347</b>	Standard head design
<b>ST3347S</b>	Standard head design with split leg
<b>ST3348</b>	Basic head design
<b>ST3348S</b>	Basic head design with split leg
<b>ST3349</b>	Plus head design
<b>ST3349S</b>	Plus head design with split leg

### Standard Features

- Safe working load and lifting capacity of 240Kg (528lbs)
- Electric height operation via integral perimeter foot switch
- Excellent height range from 46cm to 99cm
- Adjustable angle backrest from horizontal to +90°
- Adjustable head section with positive and negative angle
- Power assisted backrest and head sections incorporate anti-trap design
- Choice of standard, basic or plus head design with adjustable armrests
- Choice of standard or split leg foot section
- Breathing hole and plug fitted as standard
- Retractable wheel design with adjustable foot for stability on uneven surfaces
- Purpose made precision bearings on all pivot points
- Standard width 65cm upholstery

### Optional Accessories

- 6005** Paper Roll Holder
- 6007** Removable Head/Knee Support
- 6053** Protective Base Shroud
- 6100** Matching Operators Chair

See page 48-55 for additional information on our range of optional accessories and colour options



Model ST3347 Colour éléphant



Model ST3347S Negative head/foot section to -25°