



Width
65cm



124cm

65cm

Therapy Fixed Height Couch



The fixed height therapy couch features an extremely strong and rigid design, providing maximum patient comfort throughout all treatment with a safe working load of 240Kg (528lbs). The couch has a fully adjustable backrest, that raises from horizontal to +75° with gas strut operation. A paper roll holder, breathing hole and plug for prone treatments, are all available at an additional cost. There is an adjustable foot for added stability on one of the couch legs. The couch is supplied flat packed to reduce shipping and delivery costs.

Model

ST2700	Fixed height couch
6004	Breathing hole & plug
6005	Paper roll holder



Model ST2700 Fixed height couch with backrest partly raised

Therapy Fixed Height Couch with Drawers



A really strong and rigid couch with full length storage shelf and two storage drawers. The couch frame is powder coated for durability and has an adjustable foot on one leg for added stability. The backrest is fully adjustable by locking gas strut and a paper roll holder and breathing hole and plug are available as options. There is a wide choice of colours available, see page 55.

Models

ST2702D	Fixed height couch with storage shelf and two drawers
6004	Breathing hole & plug
6005	Paper roll holder



Model ST2702D with optional breathing hole and plug